

| BUUDY

THE COMPLETE SKINCARE GUIDE

Your Skin, *Transformed.*

A personalised guide to radiant, healthy skin using advanced LED light therapy — backed by science, designed for your life.

EXCLUSIVELY FOR BUUDY CUSTOMERS

buudy.com



Your Complete

Guide

Everything you need to know about LED light therapy, your personalised skincare routine, and how to unlock your most radiant skin.

- 01 — Welcome to the Buddy Method
- 02 — The Science of LED Light Therapy
- 03 — Your 7 Wavelengths, Explained
- 04 — Find Your Skin Type
- 05 — Your Personalised Treatment Plan
- 06 — The Daily Routine — AM & PM
- 07 — The Best Ingredients to Pair
- 08 — Complete the Ritual: Buddy IPL
- 09 — Your Results Timeline

Welcome to the

Buudy Method

You've made one of the best investments in your skin. The Buudy LED Mask isn't just a device — it's a ritual, a transformation, and a commitment to the skin you deserve.

Used by top dermatologists and aestheticians across the UK, LED light therapy is one of the most evidence-backed, non-invasive skin treatments available today. The Buudy LED Mask brings that same professional-grade technology into your home — on your schedule, at your pace.

This guide is your personal companion. We'll walk you through the science, the routines, the right products to pair, and exactly how to get the most from your device. Whether you're 25 or 65, oily or dry, acne-prone or ageing-conscious — there's a Buudy protocol for you.

—— *Your glow-up starts
here*

The Science of

LED Light Therapy

LED (Light Emitting Diode) therapy works through a biological process called photobiomodulation — the use of specific light wavelengths to trigger natural healing responses within skin cells.

When photons of light penetrate the skin, they are absorbed by chromophores — light-sensitive molecules within your cells. This triggers a cascade of benefits: increased ATP energy production, enhanced collagen synthesis, reduced inflammation, and accelerated cellular repair.

Dermis

No UV Rays

Completely safe for all skin tones and types

Non-Invasive

No pain, no downtime, no side effects


CE Certified

Approved for safe at-home use in the UK

Coll - Supproftideruis

Coll - Light

Subcutasniour Tissue



Light Meets Skin — *Cellular Magic*

1

Light Penetrates Skin

Specific wavelengths of LED light pass through the skin's surface at varying depths. Red light reaches the dermis; blue light acts on the epidermis surface.

2

Cells Absorb Energy

Mitochondria — the powerhouses of your cells — absorb light energy and become dramatically more efficient at producing ATP, the fuel for cellular repair.

3

Skin Transforms

Collagen production surges, inflammation reduces, acne bacteria are eliminated, and pigmentation fades. Results build progressively over 4–6 weeks.

Your 7 Wavelengths, *Explained*



Red

630–660 nm

Anti-ageing,
collagen boost,
firmness, reduces
fine lines



Blue

400–480 nm

Kills acne
bacteria, controls
oil, reduces
breakouts



Green

~520 nm

Pigmentation,
even skin tone,
brightens
complexion



Cyan

~490 nm

Calms
inflammation,
soothes
sensitivity,
tightens pores



Yellow

570–590 nm

Brightening,
boosts circulation,
cell renewal



Purple

Red + Blue

Dual-action: anti-
ageing and acne
simultaneously



White

400–700 nm

Full-spectrum
repair, metabolic
rejuvenation



Light

Red light at 630–660nm is absorbed by mitochondria, supercharging your cells' energy production. This stimulates fibroblasts to produce more collagen and elastin — the structural proteins that keep skin firm, plump, and youthful.

Research published in *Photomedicine and Laser Surgery* found that red LED treatment significantly reduced wrinkle depth and improved skin tone in a 12-week clinical study with UK participants.

Best For

Fine lines, wrinkles, loss of firmness, dull complexion, rosacea

Recommended Use

10 minutes, 4–5x per week. Use in the evening after cleansing.

660nm

Peak wavelength for collagen synthesis

192

High-density LEDs for full-face coverage

4–6wk

To see visible anti-ageing results

WAVELENGTH 2



Blue

Light

Blue light at 400–480nm is your skin's best weapon against acne. It penetrates the sebaceous glands and kills *Propionibacterium acnes* — the bacteria responsible for breakouts — without antibiotics or harsh chemicals.

A 2019 study in the *British Journal of Dermatology* showed a 77% reduction in inflammatory acne lesions after 8 weeks of blue LED treatment. Particularly effective for hormonal breakouts common in UK women aged 20–40.

Best For

Active breakouts, blackheads, oily skin, hormonal acne, cysts

Recommended Use

10 minutes daily during active breakouts. 3x weekly for maintenance.

460nm

Peak wavelength for killing *P. acnes* bacteria

77%

Reduction in inflammatory lesions (clinical study)

Safe to use with salicylic acid, niacinamide, and most acne topicals

Green & Yellow Light



Green — 520nm

The Brightening Wave

Green light targets melanocytes — the cells responsible for producing melanin. By regulating melanin production, it visibly fades dark spots, post-acne marks, and sun damage, while evening overall skin tone.

Particularly effective for hyperpigmentation caused by hormonal changes or sun exposure — two of the most common skin concerns for women in the UK.

Use when you have:

Dark spots, melasma, uneven skin tone, post-blemish marks, sun spots



Yellow — 570–590nm

The Glow Booster

Yellow light increases dermal blood flow and stimulates the lymphatic system, helping your skin flush toxins and deliver oxygen and nutrients more efficiently. The result? A natural, healthy flush of colour.

Ideal for the grey, tired complexion so common in the UK during winter months when natural sunlight is scarce.

Use when you have:

Dull, tired skin, poor circulation, redness, sluggish cell turnover

Cyan, Purple & White



Cyan — 490nm

The Soother

Cyan combines the calming anti-inflammatory effects of blue with green's skin-balancing qualities. It reduces redness, tightens enlarged pores, and calms reactive, sensitive skin.

Ideal for: eczema-prone, rosacea, post-procedure redness



Purple — Combined

The Dual-Action Wave

Purple combines red and blue wavelengths simultaneously — tackling ageing and acne in one treatment. Perfect for adult acne sufferers who are also concerned about fine lines and pigmentation.

Ideal for: adult acne, acne scarring, mature breakout-prone skin



White — Full Spectrum

The Regenerator

White light (400–700nm) delivers the full visible spectrum, addressing all layers of the skin simultaneously for comprehensive metabolic repair and overall rejuvenation.

Ideal for: general skin health, maintenance sessions, weekly reset

CHAPTER FOUR

Find Your

Every skin is unique. Understanding your skin type is the foundation of an effective LED routine. The right wavelength at the right time, combined with the right products, creates results that feel almost miraculous.

Normal Dry Oily Combination Sensitive Mature

*Skin
Type*

Which Skin Type Are *You*?

Dry Skin

Signs: tight feeling after cleansing, flakiness, visible lines

Dry skin often shows premature ageing as the moisture barrier is compromised. LED red light helps by boosting the production of collagen and hyaluronic acid to restore natural hydration.

Best wavelengths: Red, Yellow

Oily / Acne-Prone

Signs: shine by midday, frequent breakouts, enlarged pores

Excess sebum production feeds acne bacteria. Blue light destroys bacteria without damaging surrounding skin. Purple light adds anti-ageing benefits for those with acne scarring.

Best wavelengths: Blue, Purple, Cyan

Combination Skin

Signs: oily T-zone, dry cheeks, occasional breakouts

The most common skin type in the UK. Alternate wavelengths throughout the week — blue for oily zones, red for dry areas, and white for overall balance.

Best wavelengths: White, Blue, Red (alternate)

Sensitive Skin

Signs: redness, reactivity, rosacea, stinging from products

Start slowly with cyan and yellow to calm the skin. LED therapy is generally excellent for sensitive skin as it reduces inflammation at a cellular level without harsh chemicals.

Best wavelengths: Cyan, Yellow, Red (lower intensity)

Hyperpigmented

Signs: dark spots, melasma, uneven tone, sun damage

Green light directly targets the melanocytes responsible for dark spots. Pair with vitamin C serum post-session for amplified brightening. Yellow adds a radiance boost.

Best wavelengths: Green, Yellow, White

Mature Skin (35+)

Signs: deep lines, sagging, volume loss, crepe texture

Red light is your primary anti-ageing wavelength. Use it 5x per week. Pair with white light weekly for full-spectrum rejuvenation. Add yellow to boost circulation and radiance.

Best wavelengths: Red, White, Yellow (daily protocol)

Your Personalised *Treatment Plan*

With 7 wavelengths at your fingertips, the key to extraordinary results is knowing *when* to use each one — and how to combine them strategically throughout your week.

The protocols on the following pages are designed specifically for UK skin concerns — factoring in our climate, typical skin challenges, and lifestyle. Use them as your blueprint and adapt as your skin evolves.

10

minutes per session

3–5×

per week for best results


4–6wk

to see visible transformation

7-Day LED *Treatment Schedule*

MON	TUE	WED	THU	FRI	SAT	SUN
						
Red	Blue	Green	Red	Yellow	White	
10 min	10 min	10 min	10 min	10 min	10 min	
Anti-ageing collagen boost	Acne bacteria elimination	Fade pigmentation, brighten	Double collagen session	Pre-weekend glow boost	Full-spectrum weekly reset	
						<i>Res</i> Let your skin <i>consolidate</i> its gains

Adapt to your skin type — see personalised protocols on the following pages

A woman with dark hair tied back, wearing a white bathrobe, is looking into a large bathroom mirror. Her reflection shows her gently touching her face. The bathroom has white subway tiles. On the sink counter in the foreground, there are several skincare products, including bottles and jars. A white towel is hanging on a rack to the right.

CHAPTER SIX

The Daily *Ritual*

Great skin is built in minutes each day — not in sporadic intensive treatments. The Buddy Method creates a morning and evening ritual that takes less than 15 minutes total, but transforms your skin over time.

The order of your steps matters. Your Buddy LED Mask works best on clean, dry skin — before your serums and moisturisers. Light needs to penetrate skin directly, without barriers.

MORNING

The AM

Routine

Total: 12 minutes

1

Cleanse — 2 min

Use a gentle, pH-balanced cleanser. In the UK, hard water can strip skin — rinse with lukewarm, never hot water. Pat dry completely. Your skin must be clean and dry before LED treatment.

UK tip: Consider a filtered showerhead if you live in a hard water area (London, the South East)

2

Buddy LED Session — 10 min

Apply your Buddy LED Mask to clean, dry skin. Select your morning wavelength (see your skin type plan). Tap to start and let the light work. You can relax, listen to a podcast, or practise mindfulness.

Morning recommendation: Red (anti-ageing) or Yellow (energising glow boost)

3

Vitamin C Serum — 30 sec

Apply vitamin C serum immediately after your LED session. The photobiomodulation process increases product penetration by up to 40%, making your actives work harder. Vitamin C + LED is one of the most powerful skin brightening combinations.

4

Moisturise + SPF — 1 min

Even on overcast UK days (which is most days!), UV rays penetrate cloud cover. Apply SPF 30 or higher daily. This protects your collagen investment from UV breakdown — the single greatest cause of premature ageing.

UK essential: SPF every single day, rain or shine

EVENING

The PM

Total: 20 minutes

1

Double Cleanse — 3 min

First cleanse with an oil-based or micellar cleanser to dissolve SPF, makeup, and environmental pollution (a significant UK skin concern in cities). Second cleanse with your regular gel or foam cleanser for a truly clean base.

2

Buddy LED Session — 10 min

Your evening LED session is your most powerful treatment of the day. Skin cell regeneration peaks between 11pm–4am — your LED session primes cells for this overnight repair window.

Evening recommendation: Red (anti-ageing), Blue (acne), or Purple (dual-action)

3

Treatment Serum — 1 min

Apply your targeted treatment serum while skin is still slightly warm from LED therapy. Options: hyaluronic acid for hydration, niacinamide for pores, retinol for advanced anti-ageing. LED amplifies serum absorption significantly.

4

Night Cream / Facial Oil — 1 min

Seal everything in with a rich night cream or a few drops of facial oil. This locks in your serum actives and provides the occlusive barrier skin needs to regenerate while you sleep. For UK winters, choose a richer formula.

5

Eye Cream (Optional)

The eye area is the first to show ageing — always apply eye cream last, using your ring finger with the lightest possible pressure. Pat, don't rub. The Buddy mask covers the eye area with protective goggles.

Pro Tip from UK Dermatologists

Change your pillowcase twice weekly. Silk pillowcases reduce friction and moisture absorption — your skin retains more hydration, which amplifies your Buddy LED results overnight.

The background is a dark, textured surface. In the center, there is a large, vibrant orange slice. To its left, a large, clear water droplet is shown falling. Scattered around these elements are several molecular structures, represented by small spheres connected by lines, suggesting a scientific or chemical theme. The overall lighting is soft, highlighting the textures of the orange and the droplet.

CHAPTER SEVEN

The Best Ingredients to *Pair With LED*

LED light therapy creates a biological window of enhanced permeability immediately after each session. The right ingredients, applied within 60 seconds of removing your mask, penetrate up to 40% deeper.

Knowing which ingredients amplify your LED results — and which to avoid — is the secret most skincare guides don't tell you.

Your LED *Power Pairs*

Apply immediately after your session for maximum efficacy



Vitamin C (L-Ascorbic Acid)

Pairs with: Red, Yellow, Green

A powerful antioxidant that amplifies LED's brightening effects. Neutralises free radicals generated by UV exposure. Together with red LED, it's the gold standard for radiance. Use 10–20% concentration in the AM.



Hyaluronic Acid

Pairs with: All wavelengths

Post-LED, your skin is primed to absorb water-soluble actives. Hyaluronic acid holds 1,000x its weight in water. Apply to damp skin immediately after your session for a plumping, hydration-locking effect. A UK winter essential.



Niacinamide (Vitamin B3)

Pairs with: Blue, Cyan, Green

One of the most versatile actives in modern dermatology. Reduces pore appearance, regulates sebum, and fades pigmentation. In combination with blue LED, it's a powerhouse for acne-prone skin without irritation risk.



Retinol / Retinoids (PM Only)

Pairs with: Red (PM session only)

The gold standard anti-ageing active. Apply AFTER your evening LED session — never before, as retinol can increase photosensitivity. Red LED + retinol accelerates collagen remodelling for dramatically smoother, firmer skin over 8–12 weeks.

What to *Avoid* Before LED

Never apply these before your session — apply afterwards only

"Think of LED as an express train — the right passengers (ingredients) onboard after get first class results. The wrong ones before can derail the journey entirely."

Approach aligned with NHS dermatology recommendations

The Golden Rule

- Clean, bare skin before LED
- Serums, actives, moisturiser after LED
- Apply within 60 seconds for best absorption

YOUR DEVICE

The Buddy LED Mask — *Built Different*



192 High-Density LEDs

More than double the LED count of most competitor masks. No dead zones — every millimetre of your face and neck receives consistent, therapeutic light coverage.



Completely Cordless

2-hour rechargeable battery. No cable to manage, no bulky controller. Use it anywhere in your home — on the sofa, in bed, or even in the bath area. Freedom to relax during your treatment.



Tap Technology

Change wavelength settings with a simple tap — no fumbling with controls mid-session. Intuitive and effortless, designed to fit seamlessly into a calm, mindful routine.



Covers Neck & Jaw

Unlike most LED masks that stop at the chin, Buddy's design extends coverage to the neck and jaw — the areas that reveal age fastest. A full-coverage treatment in one session.



CE & RoHS Certified

Fully certified for use in the UK and EU. Meets all British safety standards. Dual voltage (100–240V) — use it worldwide, ideal for frequent UK travellers.

192

LEDs delivering precise wavelengths across your entire face — more coverage, more consistent results, faster transformation.

REAL RESULTS

What to Expect & *When*

WEEK 1-2

First Glow

Skin begins to look more radiant and hydrated. Pores appear slightly smaller. Active breakouts start to reduce in frequency and severity. Sleep quality often improves from the relaxation effect of sessions.

WEEK 3-4

Visible Change

Fine lines begin to soften. Dark spots start fading. Skin tone becomes noticeably more even. Others will start commenting on your skin. Collagen stimulation kicks into high gear at the cellular level.

WEEK 5-8

Transformation

The skin you've been working towards. Firmer texture, dramatically improved tone, reduced acne, and a radiant glow that holds. Many users describe feeling like they've turned back the clock by several years.

FOR BRITISH SKIN

The UK *Skin Challenges*

And exactly how Buddy addresses each one



Lack of Sunlight (Vitamin D)

British skin spends the majority of autumn and winter in almost total UV darkness. This leads to dull, grey, depleted complexions as skin cell regeneration slows without light stimulus.

Buddy solution: Yellow + White light to mimic beneficial light stimulus and boost circulation year-round



Hard Water Damage

Over 60% of England is in a hard water area. The high mineral content strips the skin's natural oils and disrupts the moisture barrier, causing dryness, sensitivity, and irritation — especially in London.

Buddy solution: Red + Cyan light to repair barrier, reduce inflammation, and restore hydration levels



Urban Pollution (London & Major Cities)

PM2.5 and nitrogen dioxide particles penetrate the skin barrier, generating free radicals that break down collagen and accelerate ageing. City skin ages faster than rural skin — often by up to 10%.

Buddy solution: Red + White light to repair cellular damage + vitamin C pairing for antioxidant protection



Hormonal Acne & Rosacea

Rosacea affects approximately 1 in 10 people in the UK. Hormonal acne spikes during stress — particularly high in British adults due to work culture and cost of living pressures. Both conditions respond excellently to LED therapy.

Buddy solution: Blue + Cyan for acne; Yellow + Cyan for rosacea; Purple for dual-concern skin

Protocol Plans by *Skin Goal*

Anti-Ageing Protocol

Mon, Wed, Fri **Red** — 10 min

Tue, Thu **Yellow** — 10 min

Saturday **White** — 10 min

Sunday **Rest**

Post-session: Retinol (PM), Vitamin C (AM)

Acne-Clearing Protocol

Mon, Wed, Fri **Blue** — 10 min

Tue, Thu **Purple** — 10 min

Saturday **Cyan** — 10 min

Sunday **Rest**

Post-session: Niacinamide, salicylic spot treatment

Brightening Protocol

Mon, Thu **Green** — 10 min

Tue, Fri **Yellow** — 10 min

Wed, Sat **White** — 10 min

Sunday **Rest**

Post-session: Vitamin C (AM), azelaic acid (PM)

Sensitive Skin Protocol

Mon, Thu **Cyan** — 10 min

Tue **Yellow** — 10 min

Saturday **Red (low)** — 8 min

Sun, Wed, Fri **Rest**

Post-session: Ceramide moisturiser, centella asiatica



"I've spent thousands on facials in London. After 6 weeks with Buddy, I cancelled my standing monthly appointment. The results are that good."

——— Sarah T., 38 —— Verified Buddy Customer, London

Make It a *Moment*

The 10 minutes of your Buddy LED session are yours. In a world of constant notifications and demands, your skincare ritual is a deliberate act of self-investment. Here's how to make the most of it — for your skin and your mind.



Put your phone down

These 10 minutes work best when you're still. The relaxation response triggered by calm + light actually enhances cellular repair. Cortisol is the enemy of healthy skin — peace is part of the protocol.



Listen to something restorative

A meditation, a sleep story, ambient music, or a relaxing podcast. Your session is a scheduled pause in the day. Protect it. Many Buddy users report improved sleep quality when they combine evening sessions with wind-down audio.



Consistency beats intensity

10 minutes, 4–5 times a week for 8 weeks will outperform sporadic marathon sessions every time. Build the habit, protect the time. Your skin will remember every session.

10

minutes a day

That's all it takes.

CHAPTER EIGHT — COMPLETE THE RITUAL

Discover the

Buudy IPL Device

You've transformed your face. Now imagine smooth, hair-free skin everywhere else — permanently, from the comfort of your home.

The Buudy IPL Hair Removal Device uses the same advanced light technology principles as your LED Mask — but instead of targeting skin cells, it targets the melanin in hair follicles to permanently reduce hair growth.



How IPL Hair Removal *Works*

IPL (Intense Pulsed Light) emits broad-spectrum light that is absorbed by the melanin pigment in your hair follicles. This light energy converts to heat, which damages the follicle's ability to regrow hair — without harming the surrounding skin.

Step 1: Flash

Broad-spectrum light flash targets the entire hair follicle and shaft. Works best on darker hair against lighter skin tones.

Step 2: Absorption

Melanin in the follicle absorbs the light energy, converting it to targeted heat that disables the follicle's growth cycle.

Step 3: Results

Over 6–8 sessions, hair becomes finer and grows back less frequently. Most users achieve 90%+ reduction in 3 months.

IPL vs Laser: What's the Difference?

Laser uses a single specific wavelength; IPL uses a broad spectrum. Both are effective for permanent hair reduction, but IPL is safer and more practical for at-home use, covering larger areas faster. Professional laser clinics in London charge £150–£300 per session. Buddy IPL delivers comparable results at home, for life.

Why *Buudy IPL* Changes Everything



Permanent Hair Reduction — At Home

The average UK woman spends £1800 per year on waxing, shaving, and threading. The Buudy IPL pays for itself in months — and delivers clinically comparable results to professional salon treatments.



Use It Anywhere

Legs, underarms, bikini line, upper lip, chin — the Buudy IPL treats every area safely. The adjustable intensity settings cater to all skin sensitivities and body areas with different hair thickness.



Also Benefits Your Skin

Beyond hair removal, IPL light pulses also stimulate collagen production in the treatment area. Many users report smoother, more even-toned skin on treated body areas as a wonderful secondary benefit.



Pain-Free Technology

Unlike waxing or traditional laser, the Buudy IPL device is designed for comfort. Adjustable intensity from 1-5 means beginners can start gently and build up. Most users describe the sensation as a mild warmth.

Your IPL *Treatment Schedule*

Follow this 8-session protocol for maximum permanent hair reduction

1

Week 1
First session

2

Week 2
Repeat

3

Week 3
Repeat

4

Week 4
Month 1 done

5

Week 6
Every 2 weeks

6

Week 8
Every 2 weeks

7

Week 10
Monthly

8

Week 14
Maintenance

Before Each Session

Shave the treatment area 24 hours before (never wax)

Cleanse skin — no body lotion or oils on area

Ensure skin is completely dry before use

Check skin tone against safety guide (Fitzpatrick scale)

After Each Session

Apply aloe vera or a fragrance-free moisturiser

Avoid sun exposure on treated areas for 24 hours

Apply SPF to any treated areas exposed to light

Avoid hot baths or saunas for 24 hours

The Complete *Buddy Routine*

Combining LED Mask + IPL Device for total body confidence



Buddy LED Mask

Face, jaw, neck

Your daily facial treatment.
Targets ageing, acne,
pigmentation, and dullness with
7 clinically proven wavelengths.
The foundation of your skin
transformation.

10 min / session — 4–5x weekly



Buddy IPL Device

Body — everywhere else

Your weekly body treatment.
Targets hair follicles for
permanent reduction on legs,
underarms, bikini, and more.
Also improves skin texture and
tone on treated areas.

Weekly sessions — 8-session initial
course



Total Confidence

Your complete transformation

Radiant, glowing face paired
with smooth, hair-free skin
everywhere else. The complete
Buddy ritual gives you
confidence in every setting —
from the office to the beach to a
black-tie event.

The full Buddy experience — face
and body

Your *Results Timeline*

LED Mask + IPL Device combined protocol



Week

1

The Beginning

Skin begins to look more radiant.
Session time becomes a cherished ritual.



Week 2-3

First Changes

Acne reduces. Skin feels smoother.

First IPL results: hair regrowth slowing, fading. Others start asking what you're



Week 4-6

Visible Glow

Tone visibly improves. Dark spots

doing.



Week 6-8

Transformation

Lines softer, texture smoother. IPL results: 50-70% hair reduction in treated areas.



Month 3

New Normal

Radiant face. Smooth body. This is your new baseline. The skin you've always wanted.



Month 6+

Lifetime Results

Maintenance sessions only. Permanent hair reduction. Anti-ageing benefits compound over time.



5 Mistakes to *Avoid*

And how to get every session working perfectly for you

0

Using it on dirty or product-covered skin

Always start with clean, dry, product-free skin. Serums, SPF, or moisturisers before your session block light penetration by up to 60%, robbing you of your results.

1

0

Skipping consistency

One long session won't undo a week of skipping. LED therapy works on cumulative exposure. 4–5 sessions of 10 minutes per week, every week. Think of it like exercise — consistency beats intensity every time.

2

0

Not pairing the right wavelength with your concern

Using blue light for anti-ageing or red light for acne won't deliver results. Match your wavelength to your concern — refer to your personalised protocol in Chapter 5. Rotate strategically through the week.

3

0

Forgetting to apply actives immediately after

The 60-second window after your session is when your skin is most receptive. Don't leave it bare. Apply your serum immediately while the biological window is open for maximum penetration.

4

0

Expecting overnight results

LED works at the cellular level — rebuilding collagen takes 28-day cell cycles. Most users see the first meaningful results at weeks 3–4. Stick with your protocol through this window. What comes next is worth the wait.

5

Frequently Asked *Questions*

Is LED light therapy safe for darker skin tones?

Yes — this is one of LED therapy's greatest advantages over lasers. LED does not target pigment and is therefore safe for all Fitzpatrick skin types (I–VI). Both red and near-infrared wavelengths carry an excellent safety profile across all ethnicities. The Buddy LED Mask is suitable for everyone in the UK, regardless of skin tone.

Can I use it during pregnancy?

LED therapy does not use UV rays or ionising radiation and is generally considered very low risk. However, as with all aesthetic devices, we recommend consulting your GP or midwife before use during pregnancy. Many healthcare professionals in the UK consider it safe, but individual guidance is always advisable.

How long until I see results?

Most users notice improved radiance and smoother texture within 1–2 weeks. More significant changes — reduced acne, softer lines, improved pigmentation — typically appear at the 4–6 week mark. Peak results come at 8–12 weeks of consistent use. Remember: LED works cumulatively, like savings in a bank. Every session counts.

Can I use it over makeup?

More Frequently Asked *Questions*

Can I use it with prescription topicals (tretinoin, antibiotics)?

LED therapy is generally compatible with most prescription skin treatments. For tretinoin or other retinoids: apply them AFTER your session in the evening only. For topical antibiotics like clindamycin: LED blue light actually enhances their efficacy. Always mention LED use to your GP or dermatologist when discussing topical prescriptions.

Does it hurt?

Not at all. LED therapy is one of the most comfortable skincare treatments available. You may feel a gentle warmth, but there is no pain, tingling, or discomfort. The Buddy mask is lightweight and contoured to fit comfortably. Many users fall asleep during their sessions — that's perfectly fine!

How does the Buddy mask compare to professional LED treatments?

Professional LED beds in UK clinics deliver higher irradiance (power per unit area) in a single session. At-home devices compensate through frequency — 4–5 sessions per week at home often produces cumulative irradiance comparable to weekly clinic visits. The Buddy mask's 192 LEDs make it one of the highest-density home devices available, closing the gap considerably.

What about IPL — is it suitable for all skin tones?

IPL works best on skin tones I–IV on the Fitzpatrick scale (light to medium-brown skin with dark hair). It is not suitable for very dark skin tones (V–VI) as the device cannot distinguish between skin and hair melanin. The Buddy IPL includes a built-in skin tone sensor for safety. Please check the Fitzpatrick scale guide included with

Seasonal Tips for *British Skin*



Autumn & Winter (Oct–Mar)

UK winters are harsh on skin: central heating dehydrates, sunlight is scarce, and cold winds damage the skin barrier. Increase your LED sessions to 5x weekly.

- Prioritise Red + White for maximum collagen and repair
- Add Yellow twice weekly to combat grey winter complexion
- Switch to a richer night cream with ceramides and fatty acids
- Layer hyaluronic acid under your serum for extra hydration



Spring & Summer (Apr–Sep)

British summers bring intermittent sunshine, occasional heatwaves, and increased UV exposure. Maintain 4x weekly sessions — adjust your routine for sun protection.

- Add Green light to address any new sun-induced pigmentation
- Wear minimum SPF 30 every single day — protect your LED investment
- Switch to lighter serums and gels in warmer months
- Blue light more frequently if summer heat triggers breakouts

Year-Round UK Essential

SPF is not optional in the UK — even in winter, UVA rays (which cause ageing) penetrate cloud cover year-round. SPF 30+ daily is the single most important step that protects and amplifies every single LED session you do.

Your Buddy Routine *By Decade*

Personalised protocols for every age and stage

20s

Prevention & clarity

Your 20s are about building habits, not fighting damage. Acne (hormonal and stress-related) is the primary concern for most UK women in this decade.

Primary: Blue (acne prevention)

Secondary: Red (early collagen building)

Weekly: White (full spectrum reset)

Post: Niacinamide, light SPF moisturiser

30s

Brightness & early prevention

The 30s bring the first fine lines, possible pigmentation from hormones or sun damage, and the start of collagen decline. This is the decade where consistency pays dividends.

Primary: Red + Green (collagen + pigmentation)

Secondary: Yellow (brightening glow)

Weekly: White (comprehensive reset)

Post: Vitamin C (AM), Retinol starter (PM)

40s

Firmness & radiance

Perimenopause begins affecting skin elasticity and hydration. Deeper lines and increased pigmentation are common. Red light becomes your primary focus, alongside intense hydration.

Primary: Red (5x weekly — non-negotiable)

Secondary: Yellow + Green (tone + glow)

Weekly: White (maximum regeneration)

Post: Retinol (PM), Peptides + HA (AM)

50s+

Restoration & confidence

Post-menopause brings significant collagen reduction — skin loses up to 30% of its collagen in the first 5 years. LED red light is the most powerful non-invasive tool available. Daily is not excessive.

Primary: Red daily (maximum anti-ageing)

Secondary: White + Yellow (full spectrum glow)

Targeted: Green (mature hyperpigmentation)

Post: Tretinoin (PM), Peptide complex + SPF50

Lifestyle Habits That *Amplify Results*

What you do between sessions matters just as much



Hydration

The UK NHS recommends 6–8 glasses of water daily. Hydrated skin is more receptive to LED therapy and heals faster. British tap water is high in minerals — consider filtered water if you live in a hard water area, as fluoride can affect cellular hydration.

Target: 2 litres daily minimum



Sleep Quality

Skin cell regeneration peaks between 11pm–4am. Getting 7–9 hours of quality sleep makes your evening LED session exponentially more effective. Growth hormone — the skin's primary repair signal — is secreted during deep sleep. Your Buddy session primes the skin; sleep completes the work.

Target: 7–9 hours before midnight



Nutrition for Skin

Collagen requires vitamin C, zinc, and amino acids to synthesise. Eat collagen-boosting foods: oily fish (omega-3s), leafy greens, citrus, eggs, berries. Limit: refined sugar (glycation damages collagen), alcohol (dehydrates and inflames), and ultra-processed foods.

Add: oily fish 2x weekly, vitamin C-rich foods daily

PROFESSIONAL PERSPECTIVE

"LED photobiomodulation represents one of the most exciting advances in non-invasive skin medicine. The evidence base for red and blue wavelengths is particularly strong. I regularly recommend at-home LED devices to my patients as a complement to — or in some cases, a replacement for — clinic treatments."

Dr. E. Richardson, Consultant Dermatologist, London

The British Association of Dermatologists recognises photobiomodulation as a clinically significant treatment modality for a range of skin conditions. LED therapy carries an excellent safety record when used as directed.

What Our Customers *Say*



"My rosacea has been completely transformed. Eight weeks in and my GP actually commented on how much better my skin looked. I've tried everything. This actually works."

Emma L., 42
Manchester — Verified Purchase



"I'm 29 and had cystic acne for 10 years. Blue light + niacinamide changed my life. I haven't had a single cyst in three months. The IPL device has also made my underarms completely smooth. Game changer."

Priya S., 29
London — Verified Purchase



"I'm 55 and post-menopausal. Red light therapy has genuinely restored firmness I thought was gone for good. I use it every night while listening to Radio 4. My ritual, my results."

Margaret H., 55
Edinburgh — Verified Purchase

The Buddy *Do's & Don'ts*



Do These Things

- Cleanse thoroughly before every session — remove all traces of makeup, SPF, and skincare
- Apply serums within 60 seconds of removing your mask to maximise the absorption window
- Wear SPF 30+ every single day — rain, shine, or overcast British skies
- Follow your personalised wavelength protocol for your specific skin concerns
- Clean your device weekly with a soft, dry cloth to maintain LED efficacy



Avoid These Mistakes

- Never use with retinol, AHAs/BHAs, or any photosensitising products applied to skin beforehand
- Don't skip sessions — inconsistency is the single biggest reason people don't see results
- Don't exceed 10–15 minutes per session — more time does not mean better results
- Don't use if you have active cold sores, open wounds, or skin infections in the treatment area
- Don't give up before 4–6 weeks — cellular collagen rebuilding takes multiple cell cycles to become visible

Build Your UK *Skincare Shelf*

The products that pair best with your Buddy LED Mask — available in the UK



Cleansers

CeraVe Hydrating Cleanser
Dry/Normal — Boots, ~£12

La Roche-Posay Effaclar Gel
Oily/Acne — Boots, ~£15

Emma Hardie Moringa Balm
Sensitive/First cleanse — £44

pH 4.5–6.5 for skin health



Serums (AM)

Timeless Vitamin C 20%
Brightening — Amazon UK, ~£20

The Ordinary HA 2% + B5
Hydration — ASOS /website, £8

Paula's Choice Niacinamide 10%
Pores/Oil — Boots, ~£38



Serums (PM)

Paula's Choice 1% Retinol
Anti-ageing — Boots, ~£54

Medik8 Crystal Retinal 3
Beginner retinoid — Medik8, £45

The Inkey List Peptide
Firming — ASOS, ~£15



SPF (Essential)

Altruist Fluid SPF 50
Best value UK SPF — Boots, £2

Ultrasun Face SPF 50+
Premium — John Lewis, ~£29

Bondi Sands Glow Serum SPF
Tinted — Superdrug, ~£18

Your Buddy on the *Go*

Don't let travel interrupt your results — tips for holidays, business trips, and city breaks

Flying Internationally

The Buddy LED Mask is dual voltage (100–240V) and works anywhere in the world with the right plug adapter. Pack it in hand luggage — it's TSA and UK security approved. A hotel session the night you arrive is the perfect jet lag antidote for your skin.

Holiday Skincare Changes

Hot, humid climates = more oil, more breakouts. Increase blue light frequency. Sunny holidays = more UV exposure. Bump your SPF to 50+ and increase green light sessions on your return to address any new pigmentation. Post-holiday red light sessions help repair any sun stress.

If You Miss Sessions

One missed week won't undo your progress. LED works cumulatively — your results bank doesn't empty overnight. Simply resume your protocol when you return. Do not try to "double up" sessions to compensate; more than 15 minutes per day is not more effective.

Travel-Sized Kit

- Buddy LED Mask + charging cable
- Universal plug adapter
- SPF 50 (100ml hand luggage)
- Travel-sized vitamin C serum
- Hyaluronic acid mist for in-flight

Worldwide ready

100–240V Dual Voltage

Works in UK, EU, USA, AUS & beyond

Buddy LED vs *The Alternatives*

Understanding the real cost and value of your skincare investment

	Buddy LED Mask	Clinic LED Bed	Monthly Facial	Prescription Topicals
Annual Cost	One-time	£1,200+/yr	£600–900/yr	£300–600/yr
Sessions per week	4–5x	1x (if budget allows)	1x monthly	Daily topicals
Wavelengths	7 wavelengths	1–3 wavelengths	Varies	N/A
Convenience	At home, anytime	Travel + appointments	Appointment needed	At home
Side effects	None	None	None–minimal	Possible irritation

Buudy IPL vs *Traditional Methods*

TRADITIONAL METHODS (PER YEAR)

Waxing (legs + underarms + bikini)	£780–£1,200
------------------------------------	-------------

Laser clinic (London pricing)	£1,500–£3,000
-------------------------------	---------------

Razors + shaving products	£100–£200
---------------------------	-----------

10-year lifetime cost	£8,000+
-----------------------	---------

BUUDY IPL DEVICE

One-time device cost	One purchase
----------------------	--------------

Sessions needed (initial course)	8 sessions
----------------------------------	------------

Maintenance sessions	Monthly (after course)
----------------------	------------------------

10-year lifetime cost	Device
-----------------------	--------

The Buudy IPL pays for itself within months of purchase — and keeps paying for life.

Stress Is the

Enemy of Your Skin

Cortisol — the body's primary stress hormone — triggers inflammation, increases sebum production, breaks down collagen, and slows wound healing. For many UK adults juggling work, family, and financial pressures, chronic stress is a constant skin saboteur.

Your Buddy LED session is a powerful antidote. 10 minutes of stillness, warm light, and intentional self-care has measurable effects on cortisol levels. Studies show that relaxation practices reduce skin inflammatory markers — your LED ritual is simultaneously a skincare and wellbeing practice.

The Buddy Mindset

This is not vanity. Investing in your skin is investing in your confidence, your wellbeing, and the way you show up in the world. You deserve to feel comfortable and confident in your own skin — and your Buddy is here to make that happen, every day.

10

minutes of intentional rest
measurably reduces cortisol

Build the habit

- Same time every day
- Pair with another ritual (tea, music)
- Track your progress in photos

Skincare *Myths, Busted*

Common misconceptions about LED therapy and skincare — set straight

MYTH

"LED therapy is just red light — they're all the same."

REALITY

Different wavelengths trigger fundamentally different cellular responses. Red builds collagen; blue destroys bacteria; green inhibits melanin; yellow boosts circulation. The 7-wavelength Buddy system gives you precision targeting that a single-wavelength device cannot match.

MYTH

"More LED = better results. Longer sessions work faster."

REALITY

LED follows a dose-response curve — there's an optimal range, and exceeding it provides no additional benefit. 10 minutes hits the therapeutic window perfectly. Frequency matters far more than session length.

MYTH

"At-home devices are just toys compared to clinic machines."

REALITY

Clinical devices have higher irradiance but are used once weekly. Home devices like Buddy are used 4–5x weekly, achieving comparable cumulative irradiance. Consistency is the deciding factor, not machine power.

MYTH

"LED doesn't work for darker skin tones."

REALITY

LED therapy is one of the few skincare technologies that works equally well across ALL skin tones. Unlike lasers or chemical peels, it doesn't target melanin. Red and near-infrared are particularly excellent and safe for Fitzpatrick types IV–VI.

A MESSAGE FROM BUUDY

Your Glow Journey

Starts Today

You've taken the first step — you invested in your skin, your confidence, and your daily ritual. Everything in this guide is designed to help you get the absolute most from that decision.

The skin of your dreams isn't a fantasy — it's the result of 10 minutes a day, the right wavelengths, the right products, and the consistency to let science do its work. You have all of that now.

Whether you're targeting acne, ageing, pigmentation, or simply want that glass-skin glow — your Buddy mask is ready. Put it on tonight. Take a photo. Then watch what happens over the next 6 weeks.

— *With love, from
Buddy.*

Ready to Add the *Buudy IPL* to Your Routine?

Your face is transforming — now imagine the same results for smooth, hair-free skin everywhere. The Buudy IPL Hair Removal Device is the perfect companion to your LED Mask routine.

- Permanent hair reduction in 8 sessions
- Works on legs, underarms, bikini, face & more
- Also improves skin texture and tone in treated areas
- Saves thousands vs. salon waxing and clinic laser

Exclusive for Buudy LED Mask Owners

Visit [buudy.com](https://www.buudy.com) to discover the Buudy IPL Device and complete your light therapy ecosystem. As a Buudy LED Mask customer, you're already part of a community committed to the best skin of their lives.





BUUDY

Your Skin Deserves

The Buddy LED Mask. The Buddy IPL Device. Your personalised routine.
This.

BUUDY.COM

This guide is for educational purposes. Always consult your dermatologist for personalised medical advice.